

User information on face masks



Warnings

- There is a risk of suffocation due to improper use.
- Do not leave children or vulnerable persons wearing the mask unsupervised.
- Its suitability during heavy physical exertion has not been tested.

Restrictions

- The face masks are intended for use exclusively during the COVID19 pandemic.
- The face masks are not a medical device.
- The face masks have not been tested for CE conformity and are intended exclusively for personal use.
- There is no guarantee or warranty for the face masks.

Application (see also: “Functional description of respiratory mask”)

- Before putting on the mask:
 - Check the face mask for a snug fit and prevent outside air from being sucked in. To increase the accuracy of fit, the mask can be adjusted in the nose area by heating. A hair dryer can be used to warm it up. Simply warm up the relevant part and then apply pressure.
 - Check for defects or damage. If the face mask is not in good condition, dispose of it and do not repair it.
 - Make sure that the valve closes correctly (for mask version with a valve).
- A face mask may only be used by one person. Do not swap the face masks with each other. It is recommended to personalize the masks by labeling them with your own name.
- The filter needs to be changed
 - if the filter becomes too moist and it's harder to breathe, or
 - after 8 hours at the latest.
- If you feel unwell while wearing a face mask, take it off immediately.
- Do not clean or reuse the filter, but always replace it with a new one.
- The selection and maintenance of the filter material is the responsibility of the user.
- Always disinfect the face mask after use.